

Running the JGI Integration Workflow

Overview

The general stages of the workflow are:

- Loading and normalizing multi-omics datasets (e.g., transcriptomics and metabolomics)
- Integrating data (quantitative values) and metadata (sample info) between datasets
- Selecting features for analysis via statistical tests
- Building and visualizing a correlation network of integrated features
- Running multi-omics factor analysis

Note: This tutorial assumes you've already completed all stages of the workflow and environment setup detailed in *setup.pdf*.

Running the Notebook

1. You should see the JupyterLab interface rendered in your browser tab.
2. Double-click the workflow notebook `/integration_workflow.ipynb` in the left menu navigator to bring it into the workspace.
3. Double-click the configuration file in `/input_data/config/project_config.yml` to bring it into the workspace. The JupyterLab interface will open the file in a text editor.
4. Run a workflow with all default parameters:
 - A. Start with the `integration_workflow.ipynb` in the workspace window.
 - B. Ensure that the kernel is "JGI Integration" by checking the top right corner of the workspace. If not, click the kernel name and select "JGI Integration" from the dropdown menu.
 - C. Run all notebook cells in order, either with the "play" button in the top menu bar or with the keyboard shortcut in each cell (ctrl/cmd/shift-click).
 - D. Each cell performs a workflow step (data loading, normalization, integration, correlation analysis, etc.) and prints some information to the standard output for review.
 - E. Review the cell print statements to see where plots and tables are saved to the `/output_data` directory; some key results or previews are also displayed in the notebook.
5. Run a workflow with customized parameters (optional):
 - A. Edit the `/input_data/config/project_config.yml` file to update workflow parameters. Use the guides provided in the `/jgi_integration/docs/*_parameters_explained.md` files to understand how parameters work.

Note: Make sure to change the `data_processing_tag` and/or `data_analysis_tag` configuration parameters – these create a new output directory to store the results and keep them distinct from previous runs. Alternatively, you can set “`overwrite=True`” in the cells that create the dataset and/or analysis object to overwrite a previous run. If you do not change the tag or overwrite, the notebook console will print an error informing you of your options.

- B. Rerun the notebook as described above from the beginning to re-load the updated configuration settings and run a full workflow with new parameters.
6. Examine the workflow results throughout the run with (A) the JupyterLab interface, (B) the notebook itself, or (C) on your local system file browser:

- A. To view the outputs within the JupyterLab web interface, navigate to the `/output_data` directory in the left side panel and click through the files and folders, or open a new workspace tab and use the built-in terminal to view files.
- B. To view the outputs within the notebook itself, you can access and view attributes of a dataset or analysis by creating a new cell and executing the command `<object>.<attribute>`. For example, running a cell with `tx_dataset.normalized_data` will show the transcriptomics count table (a dataframe) after all dataset normalization steps, or `mx_dataset.linked_metadata` will show the metabolomics metadata table after it has been linked to the other datasets.

Note: to see the names of all possible attributes that you can view for each object (dataset or analysis), run a cell with the command `vars(<object>).keys()`. For example, `vars(analysis).keys()` or `vars(mx_dataset).keys()`.

- C. To view the outputs on your system (which is mounted to the docker container), navigate to `$INTEGRATION_DIR/project_data/output_data` directory (see *setup.pdf*) and view output files as you would normally on your local filesystem.
7. Close down the docker container and JupyterLab interface with (A) the Docker Desktop app, (B) the terminal/console running the container, or (C) a new terminal/console window.
 - A. Open the Docker Desktop app and click on the Containers tab in the top left menu bar. Find the container named “jgi-integration” in the main panel and click the Stop button (black square) in the Actions column.
 - B. Open the terminal/console window that is running the container (where you ran the `... docker compose ...` command during instructions in *setup.pdf*) and press `ctrl/cmd-C` on your keyboard.
 - C. Open a new terminal/console window and stop the container with a single command:

```
> docker stop jgi-integration
```

Notes on JupyterLab

- You will not be able to delete directories on the JupyterLab file browser (the left navigation panel), so if you need to remove a directory for any reason you can do so from your local file browser (for the directories mounted to your system) or via the built-in JupyterLab terminal.
- The default working directory for JupyterLab is `/home/joyvan/work` – this is set by the docker image and should not be changed.
- Details about how to navigate the JupyterLab interface can be found on the [Jupyter documentation page](#).